



Automated SPE-method for the determination of pesticides in tea using GC-MS/MS



Lichtmanegger K., Unterluggauer H., Masselter S.


Austrian Agency for Health and Food Safety (AGES) GmbH, Competence Center for Residues of Plant Protection Products, Technikerstrasse 70, Innsbruck, Austria

Introduction

Tea crops have a worldwide economic impact and a widespread consumption due to its specific aroma as well as the promised health benefits, especially green tea. During cultivation and storage, however, various pesticides are widely used. Over the past years several reported cases of fungicide and herbicide trace level residues in green tea gained international attention and indicate the need for adequate analytical methods for routine monitoring of this commodity. The method of choice so far for monitoring purposes is the DFG-S19 multi-residue method, however with the drawback of very laborious sample preparation and strong matrix interferences. The use of tandem mass spectrometry offers various advantages in selectivity and sensitivity at low quantities and especially in such complex matrices, as it largely reduces these intrinsic matrix effects.

The proposed alternative multi-residue method presented here involves extraction with Toluene following a clean up step via siliacgel cartridges using an automatable SPE device (Liquid Handler, Gilson) coupled to GC-MS/MS analysis to provide the necessary selectivity. The scope of investigation included roughly 80 relevant pesticides (organophosphorous, organochlorine, pyrethroids, ...) and the method has been tested in several tea varieties (Green Tea, Mate Tea, Herbal tea, Black tea and in Chamomile).

Method



Results

Method validation calculations were performed on five matrix samples spiked at two concentration levels each (near the LOD of 10 µg/kg and 100 µg/kg). The results clearly demonstrate good linearity, recoveries between 70-140% for the majority of analytes and adequate precision (average RSD of 10.6%), meeting the criteria of EU guidelines (SANCO/10684/2009) (see Tab.1.).

The method was also tested for a couple of real tea samples (Tab.3) In every sample batch two spiked blank matrices (10 µg/kg /100 µg/kg) were analysed and used for calibration purposes. Only in case of MRL violations the alternative approach of standard addition for quantification has been performed.

organochlorine			pyrethroid			organophosphorous		
	RSD [%]	recovery [%]		RSD [%]	recovery [%]		RSD [%]	recovery [%]
Hexachlorbenzol	8	114	Bifenthrin	6	157	Chlorpyrifos	6	111
Dieldrin	19	123	Cypermethrin techn.	9	104	Ethion	5	111
Heptachlorepoxid	7	109	Deltamethrin	4	114	Fenitrothion	11	106
4,4'-DDD	3	109	lambda-Cyhalothrin	14	96	Malathion	10	106
2,4'-DDE	3	121	Permethrin	12	101	Pirimiphos-methyl	23	77
Endosulfan -alpha	6	115	Fenvalerate	7	121	Triazophos	8	118

Table 1: Validation data for representative Analytes in tea matrix

FAPAS (PT 2006) - tea matrix		
	assigned value (mg/kg)	analysed value (mg/kg)
p,p-DDD	0.204	0.176
Ethion	0.332	0.386
Quintozene	0.094	0.082

Table 2: Validation of the method using reference material

In addition, the presented method was successfully validated using a FAPAS tea sample from 2006 (Tab.2).

Compounds	analysed value (mg/kg)	MRL (mg/kg)
Bifenthrin	0.197	5.00
Chlorpyrifos	0.009	0.10
Cypermethrin	0.048	0.50
Fenvalerate	0.093	0.05
l-Cyhalothrin	0.07	1.00
Endosulfane	0.092	30.00

Table 3: Analysis of chinese green tea

Conclusion

The method presented demonstrates a rapid, sensitive and straightforward multi-residue approach for pesticide residue analysis in tea samples combining an automated extraction and clean-up procedure with the strength of GC-MS/MS tandem mass spectrometry and is thus suitable for food monitoring according to EU legal requirements. As a next step the method has to be validated for its suitability in daily routine analysis – with the possibility to broaden the scope of investigation.



Analytical Competence Centres